

OPSIS LiquidLINE KjelROC Scrubber

Environment-friendly Digestion



KJELROC SCRUBBER

A Scrubber is recommended for fume control and in areas where water is expensive. The KjelROC Scrubber can save costs.

RELIABLE AND EASY WITH A SMALL FOOTPRINT

The KjelROC Scrubbers are built to work in an environment with an excess of corrosive fumes and high heat. All materials are carefully selected.

- Cooler and containers in Borosilicate glass and connection tubes in Viton, in order to withstand the demanding environment and reduce maintenance costs
- Small footprint makes it easy to install the Scrubber close to the KjelROC Digestor
- Easy access to exchange containers, which are all made in high quality Borosilicate Glass.
- Accessible 6-step suction control to regulate suction power



NO WATER IS NEEDED

The KjelROC Scrubber does not require any tap water and there is a build-in pump for suction power.

OPTIONS

All KjelROC Scrubbers can be supplied with optional configuration

- Water cooler, Water Accessory to even further increase the efficiency of the fume extraction. Can also be used with activated carbon filter.
- Water Trap for digestion of liquid samples



Representantes / Distribuidores Autorizados

 Argentina

Tel: (+54 11) 5352 2500

Email: info@dastecsrl.com.ar

Web: www.dastecsrl.com.ar

 Uruguay www.dastecsrl.com.uy

 Paraguay www.dastecsrl.com.py



OPSIS AB • +46 46 72 25 00 • info@opsis.se • www.liquidline.se



LP 1018
2021 02