

GBK-1 Heat Shrink Tester



Application

To test the heating shrinkage of packaging films submerged in liquid in different temperatures. It is widely used in industries of plastic material, plastic products, petroleum, chemicals, research institute, inspection agency and so on.

Standards

ASTM D2732, GB/T 13519

Features

- ◆ Liquid tank and cover made of aluminum alloy, quick heating, warm keeping.
- ◆ LCD display temperature and timing.

Specifications

Item	Technical parameters
Sample size	≤140×140mm
Temperature range	room temperature~250℃
Temperature precision	±1℃
Instrument size	320×400×260mm
Power	1500W
Power supply	AC 220V, 50Hz

Configuration List

A、 Spare parts

NO.	Name	Remark	Quantity
1	Power cable	3*0.75 10A 1.8M	1 PC
2	Upper cover	FS-1	1 PC
3	Clamping net		1 Set
4	Clamping net bracket		2 PCS
5	Allen wrench		1 PC
6	Cross screwdriver	6-10mm	1 PC
7	Wooden box		1 PC

B、 Optional parts

NO.	Name	Unit	Quantity	Remark
1	Metrology Certificates			Optional

C、 User provide

1	Power supply	500w	220V	10A	3 holes socket	ground wire
2	Simethicone					

B、 Optional parts

NO.	Name	Unit	Quantity	Remark
1	Metrology Certificates			Optional

C、 User provide

1	Power supply	500w	220V	10A	3 holes socket	ground wire
2	Simethicone					

Note: GBPI is always committed to product innovation and improved performance, so accordingly product technical specifications are subject to change without notice. GBPI reserves the right to amend and the final power of interpretation.



Representantes / Distribuidores Exclusivos

Argentina

Tel: (+54 11) 5352 2500

Email: info@dastecsrl.com.ar

Web: www.dastecsrl.com.ar