



viscosity



Melt Flow
Indexer with
viscosity option

GÖTTFERT
 THIS IS RHEOLOGY

DASTEC

Melt Flow Indexer with viscosity option

- MFR/MVR measurement
- Viscosity option by servo motor
- Viscosity calculation via piston force
- Increased shear rate range
- Optional: Measurement of apparent flow curves according to ISO 11443
- Upper shear rate range is limited by max. force of 800 N or respectively 110 bar
- Direct operation via power panel or via PC by means of miCONNECT

